



## **GE Light Speed 16**

The LightSpeed 16 CT Scanner features an innovative solid state LightSpeed Matrix detector that enables the system to simultaneously acquire multiple adjacent slices of the anatomy. LightSpeed systems provide exceptionally clear image quality, short scan times and accurate non-invasive diagnosis for improved patient outcome. The underlying principles are the powerful 53 kW generator, the 6.3 mhu tube and the 500 millisecond scan time. In tandem with the high detection efficiency of the LightSpeed Detector, these components provide high signal output at very low patient dose.



- 36 Months \$ 10,730.30 The cost to maintain this GE Light Speed 16 is \$298.00 per day
- 60 Months \$ 8,266.26 The cost to maintain the GE Light Speed 16 is \$137.00 per day

*(The cost to maintain the GE Light Speed doesn't include electric, Ct Tech, Insurance or rent)*

**Have a question  
or need to know more?**

*Click here to contact us*



*or Call Now at  
732-262-3115*

**NJ Warehouse**  
1111 Industrial Pkwy  
Brick, NJ 08724  
Tel: (732) 262-3115

**NJ Headquarters**  
638 Mantoloking Road  
Brick, NJ 08723  
Tel: (732) 262-3115 Fax: (732) 262-3105

**Florida Office**  
11176 Creek Haven Dr.  
Riverview, FL 33569  
Tel: (732) 262-3115

**Email: [Sales@nationwideimaging.com](mailto:Sales@nationwideimaging.com)**

**[www.NationwideImaging.com](http://www.NationwideImaging.com)**